



FORMAR



PEEK CAGES		WIDTH BREITE	DEPTH TIEFE	HIGHT HÖHE	TITANIUM SIZERS
STERILE STERIL	UNSTERILE UNSTERIL				
PF-3S	PF-3	14 mm	13 mm	3 mm	SF-3
PF-4S	PF-4	14 mm	13 mm	4 mm	SF-4
PF-5S	PF-5	14 mm	13 mm	5 mm	SF-5
PF-6S	PF-6	14 mm	13 mm	6 mm	SF-6
PF-7S	PF-7	14 mm	13 mm	7 mm	SF-7
PG-3S	PG-3	16 mm	14 mm	3 mm	SG-3
PG-4S	PG-4	16 mm	14 mm	4 mm	SG-4
PG-5S	PG-5	16 mm	14 mm	5 mm	SG-5
PG-6S	PG-6	16 mm	14 mm	6 mm	SG-6
PG-7S	PG-7	16 mm	14 mm	7 mm	SG-7
PH-3S	PH-3	18 mm	14 mm	3 mm	SH-3
PH-4S	PH-4	18 mm	14 mm	4 mm	SH-4
PH-5S	PH-5	18 mm	14 mm	5 mm	SH-5
PH-6S	PH-6	18 mm	14 mm	6 mm	SH-6
PH-7S	PH-7	18 mm	14 mm	7 mm	SH-7

Titanium cages available on request / Titancages auf Anfrage erhältlich



TRAY FOR CAGES, SIZERS + TOOLS
BOX FÜR CAGES, SIZER + WERKZEUG

SMALL KLEIN	CAGE TYPES CAGE NUMMERN	PF + PG	PAZ-1
LARGE GROSS	CAGE TYPES CAGE NUMMERN	PF, PG + PH	PAZ-2



TITANIUM SETTING INSTRUMENT
TITAN SETZINSTRUMENT
180 mm
XRE-1



DISSECTING CURETTE
DISSEKTIONSKÜRETTE
6 x 200 mm
XSE-1