


FEHLING...

... THE DIFFERENCE

8

SUCTION INSTRUMENTS

SAUGINSTRUMENTE



FEHLING
INSTRUMENTS

Hanauer Landstraße 7A · 63791 Karlstein/Germany · www.fehling-instruments.de



+49 (0) 61 88 - 95 74.0



+49 (0) 61 88 - 95 74.45



info@fehling-instruments.de

SUPERLAST

MALLEABLE SUCTION CANNULAS
BIEGSAME SAUGER



Ø TUBE / ROHR	
2,0 mm	HAY-2
3,0 mm	HAY-3
4,0 mm	HAY-4


SUPERLAST instruments are malleable at regular OR temperature. They can be deformed to any helpful shape by soft bending power. While exposed to higher temperatures as i.e. that one of a washing machine or an autoclave, they automatically regain their original shape.

SUPERLAST Instrumente sind bei normaler OP-Temperatur biegsam. Sie können durch sanften Druck in die gewünschte Form gebracht werden. Unter höheren Temperaturen, wie bei der Sterillisation, nehmen sie ihre ursprüngliche Form wieder an.



MALLEABLE / BIEGBAR

Ø TUBE / ROHR	LENGTH / LÄNGE A						
	80 mm	100 mm	125 mm	145 mm	165 mm	180 mm	220 mm
0,6 mm	HAU-0	-	-	-	-	-	-
0,8 mm	HAU-4	-	-	HAI-0	-	-	-
1,0 mm	HAU-5	HAU-1	HAU-6	HAI-7	HAX-4	-	-
1,2 mm	-	HAU-2	HAU-7	HAI-8	-	-	-
1,6 mm	-	HAU-3	HAU-8	HAI-9	-	-	-
1,8 mm	-	-	-	-	HAI-1	HAI-4	HAX-1
2,0 mm	-	-	-	HAU-9	HAI-2	HAI-5	HAX-2
2,3 mm	-	-	-	-	HAI-3	HAI-6	HAX-3
3,0 mm	-	-	-	-	HAX-5	-	-



MALLEABLE / BIEGBAR

Ø TUBE / ROHR	LENGTH / LÄNGE A	
	120 mm	160 mm
1,0 mm	VOA-5	-
2,0 mm	VOA-6	VOA-1
3,0 mm	VOA-7	VOA-2
4,0 mm	VOA-8	VOA-3
5,0 mm	VOA-9	VOA-4
6,0 mm	-	VOA-0

CERAMO®

FRAZIER / FERGUSSON



CERAMO®

Ø TUBE / ROHR	TOTAL LENGTH / GESAMTLÄNGE		
	180 mm	200 mm	250 mm
1,0 mm	–	VOC-1	VOE-4
2,0 mm	VOB-2	VOC-2	VOE-5
2,5 mm	VOB-3	VOC-3	VOE-6
3,0 mm	VOB-4	VOC-4	VOE-7
4,0 mm	VOB-5	VOC-5	VOE-8
5,0 mm	VOB-6	VOC-6	VOE-9



UNCOATED
UNBESCHICHTET

Ø TUBE / ROHR	TOTAL LENGTH / GESAMTLÄNGE		
	180 mm	200 mm	250 mm
1,0 mm	–	HAR-1	–
2,0 mm	NAT-7	HAR-2	NHR-8
2,5 mm	NAT-8	HAR-0	NHS-0
3,0 mm	NAT-9	HAR-3	NHR-9
4,0 mm	NAU-1	HAR-4	NHS-1
5,0 mm	NAU-2	–	NHS-2

Ø 3,5 / 1,5 mm



JAMIESON
DISSECTION CANNULA
DISSEKTIONSSAUGER
270 mm
MMA-7



CERAMO®
SUCTION CANNULA FOR
VALVE DECALCIFICATION
ENTKALKUNGSSAUGER
240 mm
MNF-5



VALVULAR SUCTION TUBE
KLAPPENSAUGER
300 mm
MNO-0



HLM-ASPIRATOR
HLM-SAUGER
280 mm
MEL-7

FOR PEDIATRIC SURGERY
FÜR KINDERCHIRURGIE



225 mm
MCE-9



245 mm
MCE-8



210 mm
MNE-5



180 mm
SEBENING
MNE-8

DETACHABLE STRAINER
ABNEHMBARER KORB

GRECO

FOR PEDIATRIC SURGERY
FÜR KINDERCHIRURGIE



Ø TUBE / ROHR	210 mm
3,3 mm	MNL-2
5,0 mm	MNL-5



Ø TUBE / ROHR	180 mm
3,3 mm	MNL-1
5,0 mm	MNL-4



Ø TUBE / ROHR	150 mm
3,3 mm	MNL-3
5,0 mm	MNL-6

GRECO

FOR PEDIATRIC SURGERY
FÜR KINDERCHIRURGIE



$\frac{2}{3}$

DETACHABLE STRAINER / KORB ABNEHMBAR

TUBE / ROHR	180 mm
3,3 mm	MNM-1
5,0 mm	MNM-3



$\frac{2}{3}$

DETACHABLE STRAINER / KORB ABNEHMBAR

TUBE / ROHR	160 mm
3,3 mm	MNM-2
5,0 mm	MNM-4

GRECO

FOR PEDIATRIC SURGERY
FÜR KINDERCHIRURGIE



DETACHABLE STRAINER / KORB ABNEHMBAR


TUBE / ROHR	230 mm
6,0 mm	MNL-7



DETACHABLE STRAINER / KORB ABNEHMBAR


TUBE / ROHR	205 mm
6,0 mm	MNL-8

FOR PEDIATRIC SURGERY
FÜR KINDERCHIRURGIE



GRECO
HLM-CANNULA / HLM-KANÜLE

TUBE / ROHR	95 mm
3,3 mm	MNM-5
5,0 mm	MNM-7



GRECO
HLM-CANNULA / HLM-KANÜLE

TUBE / ROHR	80 mm
3,3 mm	MNM-6
5,0 mm	MNM-8



Ø 7 mm
L 25 mm
MKK-5



Ø 6 mm
L 23 mm
MKN-8

VENTRICLE SUCTION CANNULA
VENTRIKELSAUGER



¼

AORTIC DEAIRING CANNULA
AORTEN-ENTLÜFTUNGSKANÜLE

MBR-2



¼

VENTRICULAR SUCTION
CANNULA FOR 1/4" TUBE
VENTRIKELSAUGER FÜR
1/4" SCHLAUCH

MCC-0



¼

VENTRICULAR SUCTION
CANNULA
VENTRIKELSAUGER

MCC-1



Ø 2,5 mm L 48 mm
HAH-9
Ø 3,7 mm L 55 mm
HAJ-2
Ø 5,0 mm L 160 mm
HAJ-3

HEPARIN CANNULA
HEPARINKANÜLE



MALLEABLE / BIEGSAM

	Ø	1,2 mm	1,5 mm	2,0 mm
		VAD-3	VAD-5	VAD-0
		VAD-4	VAD-6	

BUTTONED CANNULAS
KNOPFKANÜLEN
80 mm